

Application Notes	This Human PRELI domain-containing protein 2 (PRELID2) ELISA Kit employs a two-site sandwich ELISA to quantitate PRELID2 in samples. An antibody specific for PRELID2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any PRELID2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PRELID2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin
-------------------	--