

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human T cell immunoreceptor with Ig and ITIM domains (TIGIT) ELISA Kit ASC-KTE62327

|                            |                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human T cell immunoreceptor with Ig and ITIM domains (TIGIT) ELISA Kit                                                                |
| Biozol Catalog Number      | ASC-KTE62327                                                                                                                          |
| Supplier Catalog Number    | KTE62327                                                                                                                              |
| Alternative Catalog Number | ASC-KTE62327-48, ASC-KTE62327-96                                                                                                      |
| Manufacturer               | Abbkine Scientific                                                                                                                    |
| Category                   | Kits/Assays                                                                                                                           |
| Application                | ELISA                                                                                                                                 |
| Species Reactivity         | Human                                                                                                                                 |
| Product Description        | This Human T cell immunoreceptor with Ig and ITIM domains (TIGIT) ELISA Kit employs a two-site sandwich ELISA to quantitate TIGIT.... |
| Range                      | Please inquire                                                                                                                        |
| Sensitivity                | Please inquire                                                                                                                        |
| Tag                        | TIGIT                                                                                                                                 |
| NCBI                       | <a href="#">201633</a>                                                                                                                |
| UniProt                    | <a href="#">Q495A1</a>                                                                                                                |
| Samples                    | Cell culture supernatantsSerumPlasmaOther biological fluids                                                                           |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>This Human T cell immunoreceptor with Ig and ITIM domains (TIGIT) ELISA Kit employs a two-site sandwich ELISA to quantitate TIGIT in samples. An antibody specific for TIGIT has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TIGIT present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TIGIT is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TIGIT bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|