

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse Thyroid Stimulating Hormone (TSH) ELISA Kit ASC-KTE70051

|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| Article Name               | Mouse Thyroid Stimulating Hormone (TSH) ELISA Kit                                                              |
| Biozol Catalog Number      | ASC-KTE70051                                                                                                   |
| Supplier Catalog Number    | KTE70051                                                                                                       |
| Alternative Catalog Number | ASC-KTE70051-48, ASC-KTE70051-96                                                                               |
| Manufacturer               | Abbkine Scientific                                                                                             |
| Category                   | Kits/Assays                                                                                                    |
| Application                | ELISA                                                                                                          |
| Species Reactivity         | Mouse                                                                                                          |
| Product Description        | This Mouse Thyroid Stimulating Hormone (TSH) ELISA Kit employs a two-site sandwich ELISA to quantitate TSH.... |
| Range                      | Please inquire                                                                                                 |
| Sensitivity                | Please inquire                                                                                                 |
| Tag                        | TSH                                                                                                            |
| Samples                    | Cell culture supernatantsSerumPlasmaOther biological fluids                                                    |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>This Mouse Thyroid Stimulating Hormone (TSH) ELISA Kit employs a two-site sandwich ELISA to quantitate TSH in samples. An antibody specific for TSH has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TSH present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TSH is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TSH bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|