

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit ASC-KTE70111

Article Name	Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit
Biozol Catalog Number	ASC-KTE70111
Supplier Catalog Number	KTE70111
Alternative Catalog Number	ASC-KTE70111-48, ASC-KTE70111-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK....
Range	Please inquire
Sensitivity	Please inquire
Tag	TNIK
NCBI	665113
UniProt	P83510
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK in samples. An antibody specific for TNIK has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNIK present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNIK is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNIK bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	--