

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ExKine(TM) Nuclear and Cytoplasmic Protein Extraction Kit ASC-KTP3001

Article Name	ExKine(TM) Nuclear and Cytoplasmic Protein Extraction Kit
Biozol Catalog Number	ASC-KTP3001
Supplier Catalog Number	KTP3001
Alternative Catalog Number	ASC-KTP3001-50, ASC-KTP3001-200
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Product Description	Abbkine ExKine(TM) Nuclear and Cytoplasmic Protein Extraction Kit enable stepwise separation and preparation of crude cytoplasmic and nuclear extracts from mammalian cultured cells or tissue....
Tag	Nuclear
Application Notes	Abbkine ExKine(TM) Nuclear and Cytoplasmic Protein Extraction Kit enable stepwise separation and preparation of crude cytoplasmic and nuclear extracts from mammalian cultured cells or tissue. This Kit is based on allowing cells to swell with hypotonic buffer. And then the cells are disrupted, the cytoplasmic fraction is removed, and the nuclear proteins are released from the nuclei by a high salt buffer. Non-denatured, active proteins are purified in less than two hours.