

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit RDD-RD-NT-PROANP-HU

| | |
|----------------------------|---|
| Article Name | Human N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit |
| Biozol Catalog Number | RDD-RD-NT-PROANP-HU |
| Supplier Catalog Number | RD-NT-ProANP-Hu |
| Alternative Catalog Number | RDD-RD-NT-PROANP-HU-96,RDD-RD-NT-PROANP-HU-48 |
| Manufacturer | reddot Biotech |
| Category | Kits/Assays |
| Application | ELISA |
| Species Reactivity | Human |
| Product Description | This Human N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit is a precoated, Sandwich ELISA assay designed to detect and quantify the target analyte in serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other... |
| Range | 39.062-2500pg/mL |
| Sensitivity | 13pg/mL |
| Samples | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids |
| Target | Human N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) |
| Application Notes | Assay Length: 4.5 hours. Experimental Method: Sandwich |