

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rabbit Myoglobin (MYO) ELISA Kit RDD-RDR-MYO-RB

|                            |                                                                                                                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Rabbit Myoglobin (MYO) ELISA Kit                                                                                                                                                                                       |
| Biozol Catalog Number      | RDD-RDR-MYO-RB                                                                                                                                                                                                         |
| Supplier Catalog Number    | RDR-MYO-Rb                                                                                                                                                                                                             |
| Alternative Catalog Number | RDD-RDR-MYO-RB-96,RDD-RDR-MYO-RB-48                                                                                                                                                                                    |
| Manufacturer               | reddot Biotech                                                                                                                                                                                                         |
| Category                   | Kits/Assays                                                                                                                                                                                                            |
| Application                | ELISA                                                                                                                                                                                                                  |
| Species Reactivity         | Rabbit                                                                                                                                                                                                                 |
| Product Description        | This Rabbit Myoglobin (MYO) ELISA Kit is a precoated, Sandwich ELISA assay designed to detect and quantify the target analyte in serum, plasma and other biological fluids . This product is for Research Use Only.... |
| Range                      | 78.1-5000pg/mL                                                                                                                                                                                                         |
| Sensitivity                | 28pg/mL                                                                                                                                                                                                                |
| Samples                    | serum, plasma and other biological fluids                                                                                                                                                                              |
| Target                     | Rabbit Myoglobin (MYO)                                                                                                                                                                                                 |
| Application Notes          | Assay Length: 3 hours. Experimental Method: Sandwich                                                                                                                                                                   |