

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse Paraoxonase 2 (PON2) ELISA Kit RDD-RDR-PON2-MU

|                            |                                                                                                                                                                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Mouse Paraoxonase 2 (PON2) ELISA Kit                                                                                                                                                                                                          |
| Biozol Catalog Number      | RDD-RDR-PON2-MU                                                                                                                                                                                                                               |
| Supplier Catalog Number    | RDR-PON2-Mu                                                                                                                                                                                                                                   |
| Alternative Catalog Number | RDD-RDR-PON2-MU-96,RDD-RDR-PON2-MU-48                                                                                                                                                                                                         |
| Manufacturer               | reddot Biotech                                                                                                                                                                                                                                |
| Category                   | Kits/Assays                                                                                                                                                                                                                                   |
| Application                | ELISA                                                                                                                                                                                                                                         |
| Species Reactivity         | Mouse                                                                                                                                                                                                                                         |
| Product Description        | This Mouse Paraoxonase 2 (PON2) ELISA Kit is a precoated, Sandwich ELISA assay designed to detect and quantify the target analyte in tissue homogenates, cell lysates and other biological fluids . This product is for Research Use Only.... |
| Range                      | 0.156-10ng/mL                                                                                                                                                                                                                                 |
| Sensitivity                | 0.052ng/mL                                                                                                                                                                                                                                    |
| Samples                    | tissue homogenates, cell lysates and other biological fluids                                                                                                                                                                                  |
| Target                     | Mouse Paraoxonase 2 (PON2)                                                                                                                                                                                                                    |
| Application Notes          | Assay Length: 3 hours. Experimental Method: Sandwich                                                                                                                                                                                          |