

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Human Serine/threonine-protein kinase receptor R3 (ACVRL1), partial (Active) BYT-ORB1882059

Article Name	Recombinant Human Serine/threonine-protein kinase receptor R3 (ACVRL1), partial (Active)
Biozol Catalog Number	BYT-ORB1882059
Supplier Catalog Number	orb1882059
Alternative Catalog Number	BYT-ORB1882059-1,BYT-ORB1882059-100,BYT-ORB1882059-20
Manufacturer	Biorbyt
Category	Proteine/Peptide
Product Description	This Recombinant Human Serine/threonine-protein kinase receptor R3 (ACVRL1), partial (Active) spans the amino acid sequence from region 22-118aa. Purity: Greater than 95% as determined by SDS-PAGE....
Molecular Weight	11.5 kDa
UniProt	P37023
Buffer	Lyophilized from a 0.2 µm filtered 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0
Source	Homo sapiens (Human)
Purity	Greater than 95% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	DPVKPSRGPLVTCTCESPHCKGPTCRGAWCTVVLVREEGRHPQEHRGCGNL HRELCRGRPTEFVNHYCCDSHLCNHNVSLVLEATQPPSEQPGTDGQ

Application Notes	<p>Biological Origin: Homo sapiens (Human). Biological Activity: Measured by its binding ability in a functional ELISA. Immobilized Human ACVRL1 at 2 µg/mL can bind Anti-ACVRL1 recombinant antibody, the EC50 is 2.417-2.971 ng/mL. Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference</p>
-------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------