

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rat Toll-like Receptor 2 (TLR-2) ELISA Kit BYT-ORB1947479

Article Name	Rat Toll-like Receptor 2 (TLR-2) ELISA Kit
Biozol Catalog Number	BYT-ORB1947479
Supplier Catalog Number	orb1947479
Alternative Catalog Number	BYT-ORB1947479-48,BYT-ORB1947479-96
Manufacturer	Biorbyt
Category	Kits/Assays
Species Reactivity	Rat
Product Description	Rat Toll-like Receptor 2 (TLR-2) ELISA Kit...
Range	0.32-20ng/mL
Sensitivity	0.19 ng/mL
UniProt	E9PTD9
Samples	serum, plasma, Tissue homogenate and Other biological samples
Target	TLR-2

Application Notes

Application Notes: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Rat TLR-2. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Rat TLR-2 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Rat TLR-2, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm. The OD value is proportional to the concentration of Rat TLR-2. You can calculate the concentration of Rat TLR-2 in the samples by comparing the OD of the samples to the standard curve