

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Phospho-Inhibitory Subunit of NF Kappa B Alpha (pIkappaBalpha) ELISA Kit BYT-ORB1948720

Article Name	Mouse Phospho-Inhibitory Subunit of NF Kappa B Alpha (pIkappaBalpha) ELISA Kit
Biozol Catalog Number	BYT-ORB1948720
Supplier Catalog Number	orb1948720
Alternative Catalog Number	BYT-ORB1948720-48,BYT-ORB1948720-96
Manufacturer	Biorbyt
Category	Kits/Assays
Species Reactivity	Mouse
Product Description	Mouse Phospho-Inhibitory Subunit of NF Kappa B Alpha (pIkappaBalpha) ELISA Kit...
Range	0.16-10ng/mL
Sensitivity	0.1 ng/mL
Samples	Tissue homogenate and Cell lysate
Target	pIkappaBalpha

Application Notes

Modification: Phosphorylated. Application Notes: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse $\text{pI}\kappa\text{B}\alpha$. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Mouse $\text{pI}\kappa\text{B}\alpha$ and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse $\text{pI}\kappa\text{B}\alpha$, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm \pm 2 nm. The OD value is proportional to the concentration of Mouse $\text{pI}\kappa\text{B}\alpha$. You can calculate the concentration of Mouse $\text{pI}\kappa\text{B}\alpha$ in the samples by comparing the OD of the samples to the standard curve