

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Legionella pneumophila subsp. pneumophila Gamma-glutamyl phosphate reductase (proA) BYT-ORB1096061

| | |
|----------------------------|--|
| Article Name | Recombinant Legionella pneumophila subsp. pneumophila Gamma-glutamyl phosphate reductase (proA) |
| Biozol Catalog Number | BYT-ORB1096061 |
| Supplier Catalog Number | orb1096061 |
| Alternative Catalog Number | BYT-ORB1096061-20,BYT-ORB1096061-100,BYT-ORB1096061-1 |
| Manufacturer | Biorbyt |
| Category | Proteine/Peptide |
| Product Description | This Recombinant Legionella pneumophila subsp. pneumophila Gamma-glutamyl phosphate reductase (proA) spans the amino acid sequence from region 1-417aa. Purity: Greater than 85% as determined by SDS-PAGE.... |
| Molecular Weight | 52.1 kDa |
| UniProt | Q5ZV33 |
| Buffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Source | Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513) |
| Purity | Greater than 85% as determined by SDS-PAGE. |
| Form | Liquid or Lyophilized powder |

| | |
|-------------------|---|
| Sequence | MNTNIANQLRAAKKATTDLNLIQSDTRTIILKTLAANLEKHIIQENQKDLSL MLEQDPYDRLLLKERILSLVNDVRKVANLPNPLGVNLFKSM PNGLSIKKIT VPLGVIAVIYESRPNVTIDIFSLCFKSGNVCILKGGKEAHFTNSYLLLLLIKNTLKNF NINTDIVCLLPPERALMTQLLNATGLVDLCIPRGSQNLINFVRDNAKIPVIETGA GIVHTYFDKSGDLGKGGKIIINNAKTRRVSCNAL |
| Application Notes | Biological Origin: Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference |