

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant *Thermus aquaticus* DNA mismatch repair protein mutL (mutL), partial BYT-ORB1096441

Article Name	Recombinant <i>Thermus aquaticus</i> DNA mismatch repair protein mutL (mutL), partial
Biozol Catalog Number	BYT-ORB1096441
Supplier Catalog Number	orb1096441
Alternative Catalog Number	BYT-ORB1096441-20,BYT-ORB1096441-100,BYT-ORB1096441-1
Manufacturer	Biorbyt
Category	Proteine/Peptide
Product Description	This Recombinant <i>Thermus aquaticus</i> DNA mismatch repair protein mutL (mutL), partial spans the amino acid sequence from region 186-312aa. Purity: Greater than 90% as determined by SDS-PAGE....
Molecular Weight	21.3 kDa
UniProt	<a href="#">P96082</a>
Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Source	<i>Thermus aquaticus</i>
Purity	Greater than 90% as determined by SDS-PAGE.
Form	Liquid or Lyophilized powder

Sequence	LFPGAGLKEAARQAFGGLLAERLFPLEKGGAFLEGLLTGPQVSRTRPDLLFLA VNGRPVALPEGVLRVRRRAYRELLPEGHYPVGVNLNLSLPPGAYRLRLDARKEE VALSKEAEAFLEEALIEEFR
Application Notes	Biological Origin: <i>Thermus aquaticus</i> . Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference