

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### SuperLight(TM) Atto HRP Reagent BAS-SLAP-200

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | SuperLight(TM) Atto HRP Reagent                                                                                                                                                                                                                               |
| Biozol Catalog Number      | BAS-SLAP-200                                                                                                                                                                                                                                                  |
| Supplier Catalog Number    | SLAP-200                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | BAS-SLAP-200                                                                                                                                                                                                                                                  |
| Manufacturer               | BioAssay Systems                                                                                                                                                                                                                                              |
| Category                   | Sonstiges                                                                                                                                                                                                                                                     |
| Product Description        | For HRP determination in ELISA, Southern Blot, Western Blot. Key Features: Single reagent: Single, stable, and ready to use. High sensitivity: Detection of HRP in the atto- to low femtogram range. Method: Luminescence, Samples: horse-radish peroxidas... |