

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Swine anti Mouse IgG (heavy and light chains), Clone: [Polyclonal], Porcine, Monoclonal NMB-SWAM/IGG(H+L)

Article Name	Swine anti Mouse IgG (heavy and light chains), Clone: [Polyclonal], Porcine, Monoclonal
Biozol Catalog Number	NMB-SWAM/IGG(H+L)
Supplier Catalog Number	SwAM/IgG(H+L)
Alternative Catalog Number	NMB-SWAM/IGG(H+L)
Manufacturer	NordicMubio
Host	Porcine
Category	Antikörper
Application	IEP, IP
Species Reactivity	Mouse
Conjugation	Unconjugated
Format	Antiserum
Target Specificity	IgG (H+L)
Cross-Adsorption (MinX)	no cross-adsorbtion
Product Description	The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in mouse serum, as...
Clonality	Monoclonal

Clone Designation	[Polyclonal]
Buffer	Neat liquid antiserum containing 0.05% sodium azide as a preservative.
Source	Purified normal IgG isolated from pooled mouse serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formula	Neat liquid antiserum containing 0.05% sodium azide as a preservative.
Antibody Type	Secondary Antibody
Application Notes	Precipitation assays. In immunoelectrophoresis use 2 μ l or equivalent against 120 μ l antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 μ l antiserum in a 3 mm diameter centre well and 2 μ l serum samples (neat and diluted) in 2 mm diameter peripheral wells.