

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Swine anti Rabbit IgG (heavy and light chains), conjugated with FITC, Clone: [Polyclonal], Porcine, Monoclonal NMB-SWAR/IGG(H+L)/FITC**

Article Name	Swine anti Rabbit IgG (heavy and light chains), conjugated with FITC, Clone: [Polyclonal], Porcine, Monoclonal
Biozol Catalog Number	NMB-SWAR/IGG(H+L)/FITC
Supplier Catalog Number	SwAR/IgG(H+L)/FITC
Alternative Catalog Number	NMB-SWAR/IGG(H+L)/FITC
Manufacturer	NordicMubio
Host	Porcine
Category	Antikörper
Species Reactivity	Mouse
Conjugation	FITC
Format	IgG
Target Specificity	IgG (H+L)
Cross-Adsorption (MinX)	no cross-adsorbtion
Product Description	The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in rabbit serum, a...
Clonality	Monoclonal
Clone Designation	[Polyclonal]

Buffer	FITC-coupled purified hyperimmune swine IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 2 ml sterile distilled water, spun dow
Source	Purified normal IgG isolated from pooled rabbit serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Form	FITC
Formula	FITC-coupled purified hyperimmune swine IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Antibody Type	Secondary Antibody
Application Notes	ELISA, Immunocytochemistry, Immunohistochemistry (frozen), (In)direct immunofluorescence.